

# FORMED ideas in motion

**NEW!**



**ARTROMOT®-K1**

CPM therapy for knee and hip joint

# ARTROMOT®-K1

CPM therapy for knee and hip joint

## Making good things great!

**ARTROMOT®-K1** combines the optimized features of the established ARTROMOT® CPM knee devices with the newest cutting-edge technology.

**ARTROMOT®-K1** has been especially developed for rental business, home care and the demanding clinical environment.

	<b>ARTROMOT®-K1</b>
Programming	Intuitive with pictograms
Stability	Metal frame architecture, maximum stability
Range of patient's size	120 – 200 cm
Weight	11 Kg
Range of motion	-10° / 0° / 120° knee 0° / 7° / 69° hip
Anatomical correct motion	Yes

## ARTROMOT®-K1 is available in different versions:

**Standard:** Offering the basics of knee CPM therapy.

**Comfort:** Offering additional special functions, including a therapy documentation and the option to combine ARTROMOT®-K1 with electrotherapy. (Ormed's ARTROSTIM®-FOCUSplus device)

**Chipcard:** Stores pre-programmed treatment parameters and the therapy documentation for each single patient.

## Indications, e.g. :

Operatively treated fractures and pseudarthrosis  
Total knee and total hip replacement • Joint mobilization under narcosis • ACL- / PCL-reconstructions • Stable osteosyntheses  
Meniscectomy • Arthroscopy • Arthrotomy

